

|                                                                                                                            |                                                                                    |                                                   |                              |
|----------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|---------------------------------------------------|------------------------------|
|  <b>DIPLOMAT<br/>DENTAL<br/>SOLUTIONS</b> | <b>Diplomat Dental s.r.o., Vrbovská cesta 17, 921 01 Piešťany, Slovak Republic</b> |                                                   |                              |
|                                                                                                                            | <b>Change description:</b>                                                         | <b>Change is valid for:</b>                       | <b>Technical information</b> |
|                                                                                                                            | <b>New Firmware Model Pro</b>                                                      | Model Pro                                         | No.: 09_26                   |
|                                                                                                                            |                                                                                    | <b>Change is valid from:</b><br><b>31/08/2026</b> | <b>Revision:</b>             |

## New firmware Model Pro

Dear valued partners, we would like to inform you about a new firmware release for our **Model Pro** line of products.

This firmware package 1.1.39, fixes the issue with **Suction control** released in **Technical information 08\_26**.

**Module CB – 1.0.0.63** – fixed a bug where upon manual connection of the unit the suction control was not working; added buffer.

Question and answers:

**Q:** I have ordered a new unit, do I need to update the firmware after installation?

**A:** No, all Model Pro units are updated to the newest available released firmware during testing.

**Q:** I have an older unit which works fine, do I need to update the unit?

**A:** If the unit works without any perceived issues then no, but as new firmware releases target specific problems, we recommend that you use the latest released firmware version.

If you have any further questions regarding this technical information or any other topic, please submit a ticket through our official [ticket.diplomat-service.com](https://ticket.diplomat-service.com) website.

Thank you,

**Your Diplomat Dental Solutions Team**