 <b>DIPLOMAT DENTAL SOLUTIONS</b>	Diplomat Dental s.r.o., Vrbovská cesta 17, 921 01 Piešťany, Slovak Republic		
	<b>Change description:</b> <b>New Diplomat Connect</b> <b>iOS app version 1.7.3</b>	<b>Change is valid for:</b> Model Pro 500/600/700/800	<b>Technical information</b> <b>No.: 024</b>
		<b>Change is valid from:</b> <b>8/1/2025</b>	<b>Revision:</b>


## New Diplomat Connect iOS app version 1.7.3

Dear valued customer, we would like to inform you of a new update to our **Diplomat Connect** app on iOS devices.

The new app version **1.7.3** has fixed the issue of previous version **1.7.0** which could have cause the Dentist profile to be save incorrectly on the server and, as a result the data was not downloaded when connecting to the unit.

The update also fixed the **Dmetec Instrument screen, Double Work Error Translation, Quick Program Selection missing** and adjusted the water bottle to the correct amount of **1300ml** instead of the **1500ml**. The update adds the new **Diplomat Dental Light** functionality and **Preset Timers**.

To automatically update the app you can turn on the auto-update function in you iOS Table in the following menu as shown in the picture.



If you have any questions regarding this technical information or any other topic, please submit a ticket through our official **[ticket.diplomat-service.com](https://ticket.diplomat-service.com)** website.

Thank you,

Diplomat Dental Solutions Team.