


| | | | |
|---|---|---|---|
|  | Diplomat Dental s.r.o., Vrbovská cesta 17, 921 01 Piešťany, Slovak Republic | | |
| | Change description: New Products and Product Configuration Updates | Change is valid for: Model Pro 500 / 600 / 700 / 800 Model One 100 / 200 / Ortho | Technical information No.: 01_26 |
| | | Change is valid from: 01/02/2026 | Revision: |

New Products and Product Configuration Updates

Dear Valued Customer, we would like to share the latest updates to our **Diplomat Dental Solutions product line**, including new **configuration options**.

We are pleased to introduce several new additions that are now available to order with our units:

New removable Cuspidor door

- Available for all carried units: Model One 100 / Model Pro 500 , 700.
- The door holds in place using metal inserts, instead of the traditional hinges. This allows for the door to be easily removed even in **tight spaces**. Clinics with limited space around the dental unit will especially benefit from this option.




New version of the legendary Woodpecker scaler

- Going forward, we are now offering the **Woodpecker V3 LED** scaler.
- This option will be available for all Model ONE / Model Pro units.

Semi-automatic disinfection of water ways for Model Pro units

- This introduces a more budget friendly option for Model Pro hygiene.

| | | | |
|---|---|---|---|
|  | Diplomat Dental s.r.o., Vrbovská cesta 17, 921 01 Piešťany, Slovak Republic | | |
| | Change description: New Products and Product Configuration Updates | Change is valid for: Model Pro 500 / 600 / 700 / 800 Model One 100 / 200 / Ortho | Technical information No.: 01_26 |
| | | Change is valid from: 01/02/2026 | Revision: |

New Integrated display option for Model Pro units

- We now offer our Model Pro line with an **integrated 9-inch capacitive touch screen** as an alternative to the tablet option.
- This makes our Model Pro units more accessible to doctors who don't have access to a tablet, or simply prefer the integrated controls more.
- The integrated option retains all functionalities of the **Diplomat Connect app**, while featuring a solid metal construction that's colour matched to the unit.



Wall charger for wireless foot controllers is now optional

- Starting from now, the wall charger for all wireless foot controllers will not be included with the units unless selected as an option in the product configuration file.
- The charger will be available to order with our units as a separate accessory in the configuration or at any time as a spare part under item number: **A653-012-901-0**

If you have any further questions regarding this technical information or any other topic, please submit a ticket through our official ticket.diplomat-service.com website.

Thank you,

Your Diplomat Dental Solutions Team